

Background

73 year old patient, male.

Non-smoker.

Pain from a mobile upper right central incisor that had received trauma approximately ten years prior to presentation.

Case performed by Dr Michael Ainsworth, UK



2. Examination revealed deep pocketing of up to 9mm mesially, grade 2 mobile, adjacent tooth grade 1 mobile and had pocket depth of 5mm distally.



4. Maryland bridge, note close adaption and support of the socket periphery.

A biphasic matrix for true bone regeneration. Built in membrane function, high graft stability.


RESORBABLE SUTURE TENTING CASE STUDY



1. Initial situation. Periapical radiograph showing resorption.



3. Atraumatic extraction leaving the socket epithelium intact.



Two weeks following extraction

 CBCT image of the extraction socket.

Ethoss Regeneration Ltd

8 Ryefield Court, Silsden, BD20 0DL, UK www.ethoss.dental • info@ethoss.dental





6. Four weeks following extraction – surgery, flap raised, completion of degranulation.

RESORBABLE SUTURE TENTING CASE STUDY



7. Palatal defect repair with EthOss®.



8. Implant placement with 2 mm cover screw.




10. Placement of 2-0 PDS tenting sutures to create "dome device" to prevent micromovement of the graft material.

A biphasic matrix for true bone regeneration. Built in membrane function, high graft stability.



9. Implant placed in ideal position for cement-retained prosthesis.




11. Placement of a "wetter" mixture of graft material in stages under the dome device, carefully drying and compressing with each stage.

Ethoss Regeneration Ltd

8 Ryefield Court, Silsden, BD20 0DL, UK www.ethoss.dental • info@ethoss.dental






12. Control radiograph at twelve weeks healing.

RESORBABLE SUTURE TENTING CASE STUDY



13. After two weeks healing post second stage.



14. Twelve months control, following crown cementation.



16. Final tissue maturation at twelve months, shows stippling, contour and volume consistent with underlying bone health.



15. Probing the distal of tooth 12 shows 1 mm pocket depth indicating adequate healing response.



17. Final tissue maturation at twelve months, shows stippling, contour and volume consistent with underlying bone health.